

ALL-OPTICAL PHYSIOLOGY WITH FAST VOLUMETRIC IMAGING IN A LIVING DROSOPHILA BRAIN

Chiao Huang,¹ Chu-Yi Tai,² Kai-Ping Yang,¹ Wei-Kun Chang,³ Kuo-Jen Hsu,^{1,3}
Yen-Yin Lin,³ Ann-Shyn Chiang,^{3,4,5,6} and Shi-Wei Chu^{1,7}

¹Department of Physics, National Taiwan University ²Institute of Biotechnology,
National Tsing Hua University ³Brain Research Center, National Tsing Hua
University ⁴Institute of Systems Neuroscience, National Tsing Hua University ⁵Institute
of Physics, Academia Sinica ⁶Kavli Institute for Brain and Mind, University of
California, San Diego ⁷Molecular Imaging Center, National Taiwan University
1, Sec 4, Roosevelt Rd, Taipei 10617, Taiwan
E-mail : swchu@phys.ntu.edu.tw

KEY WORDS: Two-photon microscopy, optogenetics, high-speed imaging, focus tunable lens, brain functional connectome

ABSTRACT

How the brain works is based on its complex neuron network. To study functional connection among neurons, compared to conventional electrophysiology, all-optical physiology [1], which uses light to stimulate and record neurons, provides attractive advantages of non-invasiveness and high spatial resolution. Although several all-optical physiology setups, including stimulation and recording techniques, have been demonstrated, the recording designs are mostly limited to single-depth imaging. Nevertheless, due to 3D distribution and fast responses of neurons, the combination of all-optical physiology and fast volumetric imaging is highly desired. In this work, simultaneous single-photon stimulation and high-speed two-photon 3D imaging is demonstrated, where the latter is achieved by integrating a tunable acoustic gradient-index (TAG) lens [2, 3]. Fully resolved functional responses are achieved within a densely distributed neuropils in a living *Drosophila* brain, as shown in Fig. 1. With precise stimulation and high spatiotemporal resolution recording, this work paves the way toward non-invasive investigation on 3D brain functional connectome *in-vivo*.




Figure 1: All-optical physiology in 3D living *Drosophila*. The functional responses in 3D coordinate are shown and activated responses are identified by generalized linear model (GLM).

REFERENCE

- [1] Emiliani, Valentina, et al. "All-optical interrogation of neural circuits." *Journal of Neuroscience*, **35**(41), 13917-13926 (2015).
- [2] Kong, Lingjie, et al. "Continuous volumetric imaging via an optical phase-locked ultrasound lens." *Nature Methods*, **12**(8), 759-762 (2015).
- [3] Hsu, Kuo-Jen, et al. "Optimizing depth-of-field extension in optical sectioning microscopy techniques using a fast focus-tunable lens." *Optics Express* **25**(14), 16783-16794 (2017).