

ZEN Intellesis: Machine Learning Based Image Segmentation

Manoj Mathew
ZEISS Research Microscopy Solutions
Carl Zeiss Pte. Ltd., 50 Kaki Bukit Place, #05-01
Singapore 415926
E-mail: manoj.mathew@zeiss.com

KEY WORDS: Machine Learning, Image Segmentation, Deep Neural Networks, Image Analysis

1. Introduction

Image segmentation is usually the biggest challenge in image analysis as segmentation lays the foundation for subsequent image analysis steps, i.e., data extraction and quantification. Besides, conventional threshold-based image-segmentation methods frequently require a specialized segmentation workflow using a combination of digital filters and tools. Moreover, there are datasets where threshold-based methods struggle or fail to generate meaningful results (e.g., images with low signal to noise ratio).

Machine learning, on the other hand, offers even non-experts the possibility to create robust and reproducible segmentation results easily. Also, machine-learning methods often yield good results with challenging images.

2. ZEN Intellesis

The software module ZEN Intellesis enables easy and precise segmentation of multidimensional images and supports images from a wide variety of different imaging modalities. ZEN Intellesis uses supervised machine learning and can be trained to automatically identify objects within an image according to a predefined set of rules (the model). In supervised learning, some features are manually labeled by the user, which forms the so-called “ground truth.” ZEN Intellesis uses pixel-based classification and offers an extremely easy to use, GUI based training interface. This enables non-experts to easily segment complex data sets without needing knowledge of machine learning. The result – a correctly classified image – can then be seamlessly integrated into the image analysis workflows to extract meaningful data. In this talk, I will describe Intellesis as a method and illustrate several examples of its applications. Figure 1 shows an example of Intellesis at work.




Figure 1: Scratch Assay of GFP expressing HeLa cells. Top left: One frame of a recorded time-lapse movie. Top right: Labeled pixels used for training of ZEN Intellesis. Bottom left: Segmentation result of the trained model – scratch area (turquoise), cell layer (orange) and mitotic cells (red). Bottom right: Based on ZEN Intellesis segmentation results, the size of the scratch area was measured over time using the ZEN Image Analysis module.