

Fluorogenic probe for fast 3D whole-cell DNA-PAINT

Kenny KH Chung^{1,2}, Zhao Zhang^{1,2}, Phylcia Kidd¹, Yongdeng Zhang¹, Nathan D Williams^{1,2}, Bennett Rollins¹, Yang Yang^{1,2}, Chenxiang Lin^{1,2}, David Baddeley^{1,2,3}, and Joerg Bewersdorf^{1,2,4}

1. Department of Cell Biology, Yale University, New Haven, CT 06520, United States

2. Nanobiology Institute, Yale University, New Haven, CT 06520, United States

3. Auckland Bioengineering Institute, University of Auckland, Auckland 1010, New Zealand

4. Department of Biomedical Engineering, Yale University, New Haven, CT 06520, United States

E-mail: kenny.chung@yale.edu

KEY WORDS: super-resolution, single-molecule switching nanoscopy, single-molecule localization microscopy.

DNA-PAINT is a single-molecule localization microscopy method based on the transient binding of DNA-based fluorescent imager probes. Unlike STORM and PALM, DNA-PAINT is resistant to photobleaching and the blinking kinetics are tunable independently of the fluorophore. However, DNA-PAINT traditionally suffers from two severe limitations: i) it is very slow, taking hours for a single image, and ii) optical sectioning is necessary to reduce background which makes it largely incompatible with standard widefield illumination.

Here we present a novel fluorogenic imager probe based on a fluorophore-quencher design. On binding to its complementary sequence, probe fluorescence increases 80 fold. To optimize for fast blinking kinetics, the probe was paired to a docking sequence with partial mismatches to decrease binding affinity. Reduced background fluorescence with our new fluorogenic probe allows imaging at significantly higher frame rates (100 Hz) compared to conventional DNA-PAINT. We demonstrate that immunolabeled microtubules in fixed COS-7 cells can be resolved as hollow tubules within minutes, with an approximately 26× faster rate of blinking. Furthermore, it is now possible to perform DNA-PAINT under widefield illumination. We demonstrated this by combining DNA-PAINT with 4Pi-SMS microscopy to image mitochondria over a thickness of more than 2 microns in 3D.

The synergy between photobleaching resistance and faster imaging rate permits the collection of super-resolution volumes at a quality and sampling-density that has not been previously achievable.

Chung et al., Fluorogenic probe for fast 3D whole-cell DNA-PAINT, bioRxiv (2020).
doi: <https://doi.org/10.1101/2020.04.29.066886>